

GraftLink<sup>™</sup> All-Inside<sup>®</sup> ACL Reconstruction with ACL TightRope<sup>®</sup> ABS

Surgical Technique



# raftLink All-Inside AC

# **GraftLink Minimally Invasive ACL Reconstruction**

The GraftLink technique provides the ultimate in anatomic, minimally invasive and reproducible ACL reconstruction.

- *Anatomic* Independent tibial and femoral socket preparation with FlipCutter<sup>®</sup> II and/or Low Profile Reamers facilitate unconstrained placement of the ACL graft.
- *Minimally Invasive* Single hamstring harvest decreases morbidity and loss of strength.\* Socket preparation with the FlipCutter II limits soft tissue dissection and preserves bone and periosteum.
- *Reproducible* The GraftLink workstation simplifies graft preparation. The tapered graft and adjustable femoral and tibial ACL TightRope buttons facilitate graft passing, fine tuning of graft depth and graft tensioning from the femoral and tibial sides.

# **Semitendinosis Harvesting and GraftLink Preparation**

In most cases only the semitendinosis is needed to create the GraftLink construct. For a less invasive option, harvest the tendon using the minimally invasive hamstring harvest technique and instruments described in technique guide LT0124.

\*Data on file





Minimally Invasive Graft Harvesting Set (AR-1279S)

# **Graft Preparation**



The GraftLink Graft Prep Attachments are placed on the base and the ACL TightRope implants are loaded into the attachments. The distance between TightRope loop ends is measured. This distance should equal 10 mm less than the desired final graft length. In this example, a TightRope is used for femoral fixation and a TightRope ABS is used for tibial fixation. Alternatively, a TightRope RT may be used for tibial fixation.



The overall graft length is measured. *Note: a length of 27 cm will yield a four-stranded GraftLink of at least 7 cm, which will provide at least 20 mm of graft in the femoral and tibial sockets.* Load the graft through the implants by folding it symmetrically over the loops. Stitch both graft ends together with a single #2 FiberLoop after passing the graft through ACL TightRopes (a). Alternatively, stitch approximately 2 cm of each graft end with one #2 FiberLoop and one #2 TigerLoop (b).



Pass one tail of each whipstitch over the graft loop and the other under the graft loop. This will ensure that the tails of the graft are tucked inside the loop during tensioning, which will facilitate tapering of ends and uniform thickness of the graft.



Once the graft is folded appropriately and the desired length is obtained, wrap the whipstitch sutures around the post to hold the construct in place.



The first stitch may now be placed. Using a "buried knot" technique, start from the inside of the graft and place the needle through the first two graft limbs.



Wrap the suture around the graft then place the needle through the second set of graft limbs from outside/in.



Tension the suture and tie a knot to secure the stitch.



This may be repeated on either end of the graft for a total of two stitches on each end.



The GraftLink Graft Prep Attachments may now be used for tensioning by simply pulling on one side until the desired tension is obtained, as read on the tensiometer. The FiberLoop whipstitch sutures may be cut off or used as supplemental fixation.



Measure the graft length and diameter. Pass both the femoral and tibial ends of the graft into the sizing block to measure diameter for socket drilling.

# **Socket Creation**

The length from the end of the femoral socket to the end of the tibial socket should be at least 10 mm longer than the graft to ensure that the graft can be tensioned fully.

Example: 70 mm graft length



Assuming a maximum intraarticular length of 30 mm, there will be approximately 20 mm of graft in the femoral and tibial socket.

Drill the femur 20 mm deep and the tibia approximately 30 mm deep to allow an extra 10 mm for tensioning.

### **Femoral Socket Preparation**

The femoral socket should be created either through the medial portal or from outside/in, using a FlipCutter II.

### Medial Portal Option



For medial portal drilling, use the TightRope Drill Pin, Transportal ACL Guides and Low Profile Reamers. Note the intraosseous length from the TightRope Drill Pin. After socket drilling, pass a suture with the TightRope Drill Pin for later graft passing.



The FlipCutter may also be used to create the femoral socket. Note the intraosseous length on the drill sleeve when pushed down to bone (a).

After "flipcutting", pass a FiberStick<sup>™</sup> suture through the Stepped Drill Sleeve and dock for later graft passing.

### **Tibial Socket Preparation**

Drill the FlipCutter into the joint. Remove the marking hook.

# **Graft Passing**



Straighten the FlipCutter blade and remove from the joint. Pass a TigerStick<sup>®</sup> into the joint and retrieve both the tibial TigerStick and the femoral FiberStick out the medial portal together with an open Suture Retriever. Retrieving both sutures at the same time will help avoid tissue interposition that can complicate graft passing. *Note: A PassPort Button Cannula™ may also be used in the medial portal to prevent tangling.* 



Flip the blade and lock into cutting position. Drill on forward, with distal traction, to cut the socket. Use the rubber ring and 5 mm markings on the FlipCutter to measure socket depth (inset).



Pass the blue button suture and the white shortening strands through the femur. Remove slack from sutures and ensure equal tension. Clamp or hold both blue and white sutures together and pull them together to advance the button out of the femur. Use markings on the loop and arthroscopic visualization of the button to confirm exit from the femoral cortex. Pull back on the graft to confirm the button is seated.



While holding slight tension on the graft, pull the shortening strands proximally, one at a time to advance the graft. Pull on each strand in 2 cm increments. *Note: The* graft can be fully seated into the femur or left partially inserted until tibial passing is complete. The latter option allows fine tuning of graft depth in each socket.



Advance the graft into the tibia by pulling on the inside of the ABS loop and whipstitch sutures.



Cinch a suture around the end of the TightRope ABS loop to use for passing (inset). Load the cinch suture and the whipstitch tails from the graft into the tibial passing suture. Pull distally on the tibial passing suture to deliver both the TightRope ABS loop and the whipstitch sutures out of the tibia distally.



Load the TightRope ABS Button onto the loop. Pull on the white shortening strands to advance the button to bone and tension graft. *Note: Ensure the button has a clear path to bone, as to not entrap soft tissue under the button.* 



Load the whipstitch sutures into the button and tie a knot for backup fixation.



# **Ordering Information**

### Implants:

ACL TightRope RT TightRope ABS	AR-1588RT AR-1588TN
TightRope ABS Button	AR-1588TB
ACL TightRope Convenience Pack	AR-1588R1S

### Instruments:

### For FlipCutter Technique:

RetroConstruction Drill Guide Set Footprint Femoral ACL Guide, right Footprint Femoral ACL Guide, left FlipCutter II, 6 mm FlipCutter II, 6.5 mm FlipCutter II, 7 mm FlipCutter II, 7 mm FlipCutter II, 8 mm FlipCutter II, 8 mm FlipCutter II, 8.5 mm FlipCutter II, 8.5 mm FlipCutter II, 9.5 mm FlipCutter II, 9.5 mm FlipCutter II, 10 mm FlipCutter II, 10 mm FlipCutter II, 11.5 mm FlipCutter II, 12 mm FlipCutter II, 13 mm	AR-1510S AR-1510FR AR-1510FL AR-1204AF-60 AR-1204AF-65 AR-1204AF-70 AR-1204AF-75 AR-1204AF-80 AR-1204AF-85 AR-1204AF-90 AR-1204AF-95 AR-1204AF-100 AR-1204AF-110 AR-1204AF-115 AR-1204AF-120 AR-1204AF-130
<u>For Medial Portal Technique:</u>	
Transportal ACL Guides (TPGs), 4 mm – 8 mm	AR-1800-04 - 08
Low Profile Reamer, 5 mm Low Profile Reamer, 5.5 mm Low Profile Reamer, 6 mm Low Profile Reamer, 6 mm Low Profile Reamer, 7 mm Low Profile Reamer, 7 mm Low Profile Reamer, 7.5 mm Low Profile Reamer, 8 mm Low Profile Reamer, 8 mm Low Profile Reamer, 8 mm Low Profile Reamer, 9 mm Low Profile Reamer, 9 mm Low Profile Reamer, 9.5 mm Low Profile Reamer, 9.5 mm Low Profile Reamer, 10 mm Low Profile Reamer, 10.5 mm Low Profile Reamer, 11.5 mm Low Profile Reamer, 11.5 mm Low Profile Reamer, 12 mm Low Profile Reamer, 12 mm Low Profile Reamer, 13 mm	AR-1405LP AR-1405LP-50 AR-1406LP-50 AR-1407LP-50 AR-1407LP-50 AR-1407LP-50 AR-1408LP-50 AR-1409LP-50 AR-1409LP-50 AR-1410LP-50 AR-1410LP-50 AR-1411LP AR-1411LP-50 AR-1412LP-50 AR-1413LP
TightRope Drill Pin, open TightRope Drill Pin, closed	AR-1595T AR-1595TC
Accessories: Suture Retriever Graft Sizing Block Graft Prep Station Base GraftLink Prep Attachment GraftLink Prep Attachment with Tensioner Suture Cutter for ACL TightRope	AR-12540 AR-1886 AR-2950 AR-2951-1 AR-2951-2 AR-4520
Suture:	
0 FiberWire, 38" (blue) w/Tapered Needle, 22.2 mm 1/2 circle FiberStick, #2 FiberWire, 50" (blue) one end stiffened	AR-7250 AR-7209
TigerStick, #2 TigerWire, 50" (white/black) one end stiffened #2 FiberLoop w/Straight Needle, 20" (blue)	AR-7209T
76 mm needle w/7 mm loop	AR-7234
#2 HgerLoop w/Straight Needle, 20" w/ HgerWire (white/green), 76 mm needle w/7 mm loop	AR-7234T

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use.



U.S. PATENT NOS. 6,716,234; 7,029,490; 7,147,651 and PATENT PENDING

©2012, Arthrex Inc. All rights reserved. LT1-0157-EN\_C